








THORATEC
CORPORATION

FOR USE ANYTIME

CONTACT: Susan Benton Russell
Phone: 310.697.3488
susan@bentoncommunications.com

Thoratec Mechanical Circulatory Support (MCS) Product Portfolio

					
Approved Indications (U.S.)	CentriMag [®] + Blood Pump	Thoratec [®] PVAD™	Thoratec [®] IVAD™	HeartMate [®] XVE	HeartMate II [®]
Bridge-to-Transplantation (BTT)		✓	✓	✓	✓
Post-Cardiotomy Recovery (PC)		✓	✓		
Home Discharge		✓	✓	✓	✓
Destination Therapy (DT)				✓	✓
Extracorporeal Circulatory Support	✓	paracorporeal			
Planned Indications (U.S.)					
Post-Cardiotomy Recovery (PC)	✓				
Product Attributes					
Acute Support (hours - weeks)	✓ (up to 6 hours)	✓			
Intermediate Support (months)	✓ RVAD use [†]	✓	✓	✓	✓
Chronic Support (years)			✓	✓	✓
Bilateral (right & left heart support)	✓	✓	✓		
Mobility		✓	✓	✓	✓
Small Patient Size	✓	✓	✓		✓
Pneumatic		✓	✓		
Electric	✓			✓	✓
Pulsatile		✓	✓	✓	
Rotary					✓
Centrifugal	✓				

B100-1109

*CentriMag is a registered trademark of and is manufactured by Levitronix, LLC.

† Humanitarian Use Device. The Levitronix CentriMag RVAS is authorized by Federal law to provide temporary circulatory support for up to 30 days for patients in cardiogenic shock due to acute right ventricular failure. The effectiveness of this device for this use has not been demonstrated. Distribution of this device is restricted to use by or on the order of a physician.